



CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
www.lab-honeywell.com

Brand: Honeywell Riedel-de Haën
Product: 34858
Lot No: J162A
Production date: 11.Jun.2019
Best before: 05.Mar.2034

Ethyl acetate CHROMASOLV, for HPLC, 99.7%

Parameter	Specification	Units	Result	Units
assay (GC)	min. 99.7	%	99.99	%
non-volatile matter	max. 0.0005	%	0.00005	%
water (Karl Fischer)	max. 0.03	%	0.002	%
free acid (as CH₃COOH)	max. 0.005	%	<0.005	%
absorbance at 400 nm	max. 0.01		<0.001	
absorbance at 300 nm	max. 0.01		0.001	
absorbance at 270 nm	max. 0.10		0.004	
absorbance at 260 nm	max. 0.30		0.095	
absorbance at 255 nm	max. 0.70		0.530	



QC Release Date: 19.Jun.2019

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell

Quality Management

This document was produced electronically and is valid without a signature