



CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
lab.honeywell.com

Brand: Honeywell Riedel-de Haën
Product: 34480
Lot No: O214A
Production date: 01.Aug.2024
Expiry Date: 22.Jul.2026

Acetone CHROMASOLV TM, for pesticide residue analysis

Parameter	Specification	Units	Result	Units
assay (GC)	min. 99.8	%	99.98	%
non-volatile matter	max. 0.0005	%	<0.0005	%
water (Karl Fischer)	max. 0.2	%	0.03	%
pesticides filtered via 0,2µm pore size filter	max. 5	ng/l	<5	ng/l

QC Release Date: 08.Aug.2024

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell
Quality Management

This document was produced electronically and is valid without a signature