

# Honeywell

## CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH  
 Wunstorferstrasse 40  
 Seelze, 30926  
 Telefon: +49 5137 999-0  
[lab.honeywell.com](http://lab.honeywell.com)

Brand: Honeywell Fluka  
 Product: 03124  
 Lot No: O1490  
 Production date: 28.May.2024  
 Rec. Retest Date: 18.May.2026

Potassium iodide puriss., meets analytical specification of Ph. Eur., BP, USP

Parameter	Specification	Units	Result	Units
assay (calc. to the dried substance)	99.0-100.5	%	99.4	%
loss on drying (105°C)	max. 0.5	%	0.03	%
alkalinely reacting substances	complying		complying	
appearance of the solution	complying		complying	
iron (Fe)	max. 10	ppm	<10	ppm
heavy metals (as Pb)	max. 10	ppm	<10	ppm
iodate (JO3)	complying		complying	
iodate (IO3)	max. 0.0004	%	<0.0004	%
sulfate (SO4)	max. 100	ppm	<100	ppm
thiosulfate (S2O3)	max. 0.01	%	<0.01	%
nitrate, nitrite, ammonium	complying		complying	
thiosulfate (S2O3) und Ba	complying		complying	
residual solvents	complying		complying	



QC Release Date: 06.Jun.2024

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell

Quality Management

This document was produced electronically and is valid without a signature